acumed®

Anatomic Radial Head System

MediMark® Europe Sarl.  
11 rue Emile ZOLA. BP 2332  
38033 GRENOBLE CEDEX 2  
FRANCE  
+33 4 76 86 43 22

Acumed® LLC  
5885 NW Cornelius Pass Road  
Hillsboro, OR 97124-9432  
+1 (503) 627-9957  
www.acumed.net

PKGI-42-E  
EFFECTIVE 4/2015

STERILE R

REF

LOT

CE

0473

\* P K G I - 4 2 - E \*

OBSOLETE  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

# TABLE OF CONTENTS

Click on language to navigate to page.

English – US .....	3
Dansk – DA .....	9
English – EN .....	15
Deutsch – DE .....	21
Ελληνικά – EL .....	28
Español – ES .....	35
Français – FR .....	42
Italiano – IT .....	49
Nederlands – NL .....	56
Norsk – NO .....	63
Português – PT .....	69
Suomi – FI .....	76
Svenska – SV .....	83
Türkçe – TR .....	89

**OBSOLETE**  
Visit [www.acumed.net/ifu](http://www.acumed.net/ifu)  
for the latest version.

# US ACUMED® ANATOMIC RADIAL HEAD SYSTEM

## FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

**DESCRIPTION:** The Acumed Anatomic Radial Head System of implants and instruments are designed to replace the articular surface of the radial head.

**INFORMATION FOR USE:** Physiological dimensions limit the sizes of implant appliances. The surgeon must select the type and size that best meets the patient's requirements for close adaptation and firm seating with adequate support.

Many of Acumed Instruments and Accessories are reusable; however they have a limited life span. Prior to and after each use, the instruments must be inspected where applicable for sharpness, wear, damage, proper cleaning, corrosion and integrity of the connecting mechanisms. Particular care should be paid to hex drivers, drill bits and instruments used for cutting or implant insertion.

Although the physician is the learned intermediary between the company

and the patient, the important medical information given in this document should be conveyed to the patient.

**SURGICAL TECHNIQUES:** Surgical techniques are available describing the uses of this system. It is the responsibility of the surgeon to be familiar with the procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use.

**INDICATIONS:** The Acumed Anatomic Radial Head System and accessories are designed specifically for (1) Replacement of the radial head for degenerative or post-traumatic disabilities presenting pain, crepitation, and decreased motion at the radio-humeral and/or proximal radio-ulnar joint with: joint destruction and/or subluxation, resistance to conservative treatment, (2) Primary replacement after fracture of the radial head, (3) Symptomatic sequelae after radial head resection,

(4) Revision following failed radial head arthroplasty.

**CONTRAINDICATIONS:** Contraindications for the system are active or latent infection; sepsis; osteoporosis, insufficient quantity or quality of bone or soft tissue and material sensitivity. If material sensitivity is suspected, tests should be performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices. These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

**WARNINGS:** For safe effective use of the implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for the device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading

associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical

evaluation using MRI equipment<sup>1</sup>.

For safe effective use of any Acumed instrument, the surgeon must be familiar with the instrument, the method of application, and the recommended surgical technique. Instrument breakage or damage can occur when an instrument is subjected to excessive loads, speeds, or dense bone. The patient must be cautioned, preferably in writing as to the risks associated with these types of instruments.

<sup>1</sup> Shellock, F. G. Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition, Biomedical Research Publishing Group, 2011.

**PRECAUTIONS:** An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking. Such stress concentrations can lead to failure.

Mixing implant components from different manufacturers is not recommended for metallurgical, mechanical and functional reasons.

The benefits from implant surgery may not meet the patient's expectations or may

deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

**ADVERSE EFFECTS:** Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur. Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

#### **CLEANING:**

**Implant Cleaning:** This product is provided sterile, and should not be re-cleaned. Acumed does not recommend re-cleaning or re-sterilization of sterile-

packaged product.

**Instrument Cleaning:** Acumed Instruments and Accessories must be thoroughly cleaned before reuse, following the guidelines below:

**Warnings & Precautions**

- Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Do not allow contaminated instruments to dry prior to cleaning/reprocessing. Excess blood or debris should be wiped off to prevent it from drying onto the surface.
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.
- Do not use metal brushes or scouring pads during manual cleaning process.
- Use cleaning agents with low foaming surfactants for manual cleaning in order to see instruments in the cleaning solution. Cleaning agents must be easily rinsed from instruments to prevent residue.
- Mineral oil or silicone lubricants should

not be used on Acumed instruments.

- Neutral pH enzymatic and cleaning agents are recommended for cleaning reusable instruments. It is very important that alkaline cleaning agents are thoroughly neutralized and rinsed from instruments.
- Surgical instruments must be dried thoroughly to prevent rust formation, even if manufactured from high grade stainless steel.
- All instruments must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.
- Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.

**Manual Cleaning/Disinfection Instructions**

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer.

Fresh solutions should be prepared when existing solutions become grossly contaminated.

2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of twenty (20) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush. For exposed springs, coils, or flexible features: Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.
3. Remove the instruments and rinse thoroughly under running water for a minimum three (3) minutes. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.

4. Place the instruments, fully submerged, in an ultrasonic unit with cleaning solution. Actuate all moveable parts to allow detergent to contact all surfaces. Sonicate the instruments for a minimum of ten (10) minutes.
5. Remove the instruments and rinse in deionized water for a minimum of three (3) minutes or until all signs of blood or soil are absent in the rinse stream. Pay special attention to cannulations, and use a syringe to flush any hard to reach areas.
6. Inspect instruments under normal lighting for the removal of visible soil.
7. If visible soil is seen, repeat the sonication and rinse steps above.
8. Remove excess moisture from the instruments with a clean, absorbent, nonshedding wipe.

**Combination Manual/Automated Cleaning and Disinfecting Instructions**

1. Prepare enzymatic and cleaning agents at the use-dilution and temperature recommended by the manufacturer. Fresh solutions should be prepared when existing solutions become grossly contaminated.

2. Place instruments in enzymatic solution until completely submerged. Actuate all moveable parts to allow detergent to contact all surfaces. Soak for a minimum of ten (10) minutes. Use a nylon soft bristled brush to gently scrub instruments until all visible debris is removed. Pay special attention to hard to reach areas. Pay special attention to any cannulated instruments and clean with an appropriate bottle brush.

*Note: Use of a sonicator will aid in thorough cleaning of instruments. Using a syringe or water jet will improve flushing of difficult to reach areas and any closely mated surface.*

3. Remove instruments from enzyme solution and rinse in deionized water for a minimum of one (1) minute.
4. Place instruments in a suitable washer/disinfecter basket and process through a standard washer/disinfecter cycle. The following minimum parameters are essential for thorough cleaning and disinfection.

Step	Description
1	Two (2) minute prewash with cold tap water
2	Twenty (20) second enzyme spray with hot tap water
3	One (1) minute enzyme soak
4	Fifteen (15) second cold tap water rinse (X2)
5	Two (2) minute detergent wash with hot tap water (64-66°C/146-150°F)
6	Fifteen (15) second hot tap water rinse
7	Ten (10) second purified water rinse with optional lubricant (64-66°C/146-150°F)
8	Seven (7) minute hot air dry (116°C/240°F)
<i>Note: Follow washer/disinfecter manufacturer's instructions explicitly</i>	

**Automated Cleaning/Disinfection Instructions**

- Automated washer/dryer systems are not recommended as the only cleaning method for surgical instruments.
- An automated system may be used as a follow up process to manual cleaning.
- Instruments should be thoroughly inspected prior to sterilization to ensure effective cleaning.

**STERILITY:**

System components may be provided sterile or nonsterile.

**Sterile Product:** Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilization of a device provided sterile should be handled according to hospital protocol. Acumed does not recommend resterilization of sterile-packaged product.

**Non-Sterile Product:** Unless clearly labeled as sterile and provided in an unopened sterile package provided by Acumed, all implants and instruments must be considered nonsterile, and sterilized by the hospital prior to use. Nonsterile devices have been validated using the sterilization parameters listed below, in fully-loaded trays with all parts placed appropriately.


**Sterilization Methods**

- Consult your equipment manufacturer’s written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN “Recommended Practices for Sterilization in Perioperative Practice Settings” and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed

according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

**STORAGE INSTRUCTIONS:** Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

<b>80-2000/80-2001 These trays have been validated to the parameters below:</b>	
<b>Gravity Displacement Autoclave: NOT RECOMMENDED</b>	
<b>Pre-Vacuum Autoclave:</b>	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Dry Time:	30 minutes
<b>TR-0002/TR-0003 This tray has been validated to the parameters below:</b>	
<b>Gravity Displacement Autoclave:</b>	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	60 minutes
Dry Time:	125 minutes
<b>Pre-Vacuum Autoclave:</b>	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	25 minutes
Dry Time:	70 minutes

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Catalogue number
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

**APPLICABILITY:** These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

**FURTHER INFORMATION:** To request further material, please see the contact information listed on this document.

Visit [www.acumed.net/ifu](http://www.acumed.net/ifu) for the latest version.

**Cautions: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital. Professional Use Only.**

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assurance in health care facilities.

- Lyn-sterilisering anbefales ikke, men hvis den anvendes, bør den kun udføres i følge bestemmelserne i ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

<b>80-2000/80-2001 Disse bakker er blevet godkendt iht. parametrene angivet herunder:</b>	
<b>Autoklaving, tyngdekraftsforskydning: FRARÅDES</b>	
<b>Autoklaving, prævakuum:</b>	
Tilstand:	Indpakket
Eksponeringstemperatur:	132° C (270° F)
Eksponeringsperiode:	4 minutter
Tørreperiode:	30 minutter
<b>TR-0002/TR-0003 Disse bakker er blevet godkendt iht. parametrene angivet herunder:</b>	
<b>Autoklaving, tyngdekraftsforskydning:</b>	
Tilstand:	Indpakket
Eksponeringstemperatur:	132° C (270° F)
Eksponeringsperiode:	60 minutter
Tørreperiode:	125 minutter
<b>Autoklaving, prævakuum:</b>	
Tilstand:	Indpakket
Eksponeringstemperatur:	132° C (270° F)
Eksponeringsperiode:	25 minutter
Tørreperiode:	70 minutter

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associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/ or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/ or treatment to fail. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/ or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical

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**ADVERSE EFFECTS:** Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/ or loosening may occur. Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

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Step	Description
1	Two (2) minute prewash with cold tap water
2	Twenty (20) second enzyme spray with hot tap water
3	One (1) minute enzyme soak
4	Fifteen (15) second cold tap water rinse (X2)
5	Two (2) minute detergent wash with hot tap water (64-66°C/146-150°F)
6	Fifteen (15) second hot tap water rinse
7	Ten (10) second purified water rinse with optional lubricant (64-66°C/146-150°F)
8	Seven (7) minute hot air dry (116°C/240°F)

*Note: Follow washer/disinfecter manufacturer's instructions explicitly*

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
#### **Sterilization Methods**

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed

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<b>80-2000/80-2001</b> These trays have been validated to the parameters below:	
<b>Gravity Displacement Autoclave:</b> NOT RECOMMENDED	
<b>Pre-Vacuum Autoclave:</b>	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Dry Time:	30 minutes
<b>TR-0002/TR-0003</b> This tray has been validated to the parameters below:	
<b>Gravity Displacement Autoclave:</b>	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	60 minutes
Dry Time:	125 minutes
<b>Pre-Vacuum Autoclave:</b>	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	25 minutes
Dry Time:	70 minutes

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