ACUMED®

DUAL-TRAK™ FIXATION SCREW SYSTEMS

MediMark® Europe Sarl.
11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

Acumed®
5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1 (503) 627-9957
www.acumed.net

PKGI-62-A
Effective 5/2012

STERILE R

REF

LOT

0473

* P K G I - 6 2 - A *

DUAL-TRAK™ FIKSERINGSKRUESYSTEMER

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WARNINGS: For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, non-union, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant including the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, non-union, or incomplete healing including possible nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The implant system has not been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment, nor has it been tested for heating or migration in the MR environment.

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ADVERSE EFFECTS: Fracture of the implant due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion, implant migration and/or loosening. Metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material. Pain, discomfort, or abnormal sensations due to the presence of an implant. Nerve or soft tissue damage may result from the presence of an implant or to surgical trauma. Necrosis of bone or bone resorption. Necrosis of the tissue or inadequate healing.

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Gravity Displacement Autoclave	Full Cycle Time: 270° F (132°C) for 60 minutes Dry Time: 125 minutes
Pre-Vacuum Autoclave	Full Cycle Time: 270° F (132° C) for 25 minutes Dry Time: 70 minutes

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- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

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DUAL-TRAK™ FIXATION SCREW SYSTEMS

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: The Acumed fixation devices are designed to provide fixation of various fractures and osteotomies while they heal.

INFORMATION FOR USE: The surgeon must select the type and size implant that best meets the patient's surgical needs.

INDICATIONS: The Dual-Trak™ Clavicle Screw is intended for fusion and fixation of the clavicle.

CONTRAINDICATIONS: Active or latent infection. Osteoporosis, insufficient quantity or quality of bone/soft tissue. Material sensitivity. If suspected, tests are performed prior to implantation. Sepsis. Patients who are unwilling or incapable of following postoperative care instructions. This device is not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

WARNINGS: For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device. The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant including the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing including possible nerve or soft tissue damage related to either surgical trauma or the presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The implant system has not been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment, nor has it been tested for heating or migration in the MR environment.

PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking. Such stress concentrations can lead to failure.

ADVERSE EFFECTS: Fracture of the implant due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion; implant migration and/or loosening. Metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material. Pain, discomfort, or abnormal sensations due to the presence of an implant. Nerve or soft tissue damage may result from the presence of an implant or to surgical trauma. Necrosis of bone or bone resorption. Necrosis of the tissue or inadequate healing.

STERILITY: This product is provided sterile and non-sterile. Sterile device was exposed to a minimum dose of 25.0-kGy gamma irradiation. Sterilization may be performed using one of the following methods.

STERILIZATION METHODS

FOR ACUMED UNIVERSAL TRAY WITH DUAL-TRAK™ PLATTER (80-0341 WITH 80-0785):

Gravity Displacement Autoclave	Full Cycle Time: 270° F (132° C) for 60 minutes Dry Time: 125 minutes
Pre-Vacuum Autoclave	Full Cycle Time: 270° F (132° C) for 25 minutes Dry Time: 70 minutes

FOR DUAL-TRAK™ PLATTER & LID ONLY (80-0785 WITH 80-0786):

Gravity Displacement Autoclave	NOT RECOMMENDED
Pre-Vacuum Autoclave	Full Cycle Time: 270° F (132° C) for 4 minutes Dry Time: 30 minutes

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN® Recommended Practices for Sterilization in Perioperative Practice Settings® and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

CAUTION: For Professional Use Only.

CAUTION: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital.