ACUMED®

ACUMED® PELVIC BONE PLATE SYSTEM

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

MediMark® Europe Sarl.

11 rue Emile ZOLA. BP 2332
38033 GRENOBLE CEDEX 2
FRANCE
+33 4 76 86 43 22

EC REP

Acumed®

5885 NW Cornelius Pass Road
Hillsboro, OR 97124-9432
+1 (503) 627-9957
www.acumed.net

PKGI-63-B
EFFECTIVE 7/2013

STERILE R

REF

LOT

0473

* P K G I - 6 3 - B *

TABLE OF CONTENTS

Click on language to navigate to page.

English – US.....	3
Dansk – DA.....	8
English – EN.....	13
Deutsch – DE.....	18
Ελληνικά – EL.....	24
Español – ES.....	30
Français – FR.....	35
Italiano – IT.....	41
Nederlands – NL.....	47
Norsk – NO.....	53
Português – PT.....	58
Suomi – FI.....	65
Svenska – SV.....	68
Türkçe – TR.....	73

OBSOLETE
Visit www.acumed.net/ifu
for the latest version.

US ACUMED® PELVIC BONE PLATE SYSTEM

FOR THE PERSONAL ATTENTION OF THE OPERATING SURGEON

DESCRIPTION: The Acumed Pelvic Bone Plate System of bone plates, screws and accessories is designed to provide fixation for fractures, fusions or osteotomies.

The Acumed Pelvic Bone Plate System also includes instruments to facilitate placement of implants.

INFORMATION FOR USE: Physiological dimensions limit the sizes of implant appliances. The surgeon must select the type and size that best meets the patient's requirements for close adaptation and firm seating with adequate support.

Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

SURGICAL TECHNIQUES: Surgical techniques are available describing the uses of this system. It is the responsibility of the surgeon to be familiar with the

procedure before use of these products. In addition, it is the responsibility of the surgeon to be familiar with relevant publications and consult with experienced associates regarding the procedure before use.

INDICATIONS: The Acumed Pelvic Bone Plate system includes plates, screws and accessories designed to provide fixation during fractures, fusions, and osteotomies for the acetabulum, sacrum, ilium, and entire pelvic ring, as well as treatment of sacroiliac joint dislocations and symphysis pubis disruptions.

CONTRAINDICATIONS: Contraindications for the system are active or latent infection, sepsis, osteoporosis, insufficient quantity or quality of bone/soft tissue, and material sensitivity. If sensitivity is suspected, tests are performed prior to implantation. Patients who are unwilling or incapable of following postoperative care instructions are contraindicated for these devices. These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

WARNINGS: For safe effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique for this device.

The device is not designed to withstand the stress of weight bearing, load bearing, or excessive activity. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening or migration. The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant including the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or weight bearing or load bearing, particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete healing including possible nerve or soft tissue damage related to either surgical trauma or the

presence of the implant. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail. The components of these systems have not been tested for safety, heating, or migration in the MRI environment. The implants may cause distortion and/or block the view of anatomic structures on radiographic images. Similar products have been tested and described in terms of how they may be safely used in post-operative clinical evaluation using MRI equipment¹.

1 Sherlock, F. G. Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition. Biomedical Research Publishing Group, 2011.

PRECAUTIONS: An implant shall never be reused. Previous stresses may have created imperfections, which can lead to a device failure. Instruments shall be inspected for wear or damage prior to usage. Protect implants against scratching and nicking. Such stress concentrations can lead to failure.

Mixing implant components from different manufacturers is not recommended for metalurgical, mechanical and functional reasons.

The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

ADVERSE EFFECTS: Possible adverse effects are pain, discomfort, or abnormal sensations and nerve or soft tissue damage due to the presence of an implant or due to surgical trauma. Fracture of the implant may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur. Metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material may occur. Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing may result from the presence of an implant or due to surgical trauma.

CLEANING:

Implant Cleaning: Implants should not be re-used. However, any implant that has not been used, but has

become soiled, should be processed according to the following:

Warnings & Precautions

- Any implant contaminated with blood, tissue and/or bodily fluids/matter should be processed according to hospital protocol.
- Do not use an implant if the surface has been damaged. Damaged implants should be discarded.
- Users should wear appropriate personal protective equipment (PPE)
- All users should be qualified personnel with documented evidence of training and competency. Training should be inclusive of current applicable guidelines and standards and hospital policies.

Manual Processing

Equipment: Soft bristled brush, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5.

1. Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and

concentration.

- Carefully wash the implant manually. Do not use steel wool or abrasive cleaners on implants.
- Rinse implant thoroughly with DI or purified water. Use DI or purified water for final rinse.
- Dry the implant using a clean soft, lint-free cloth to avoid scratching the surface.

Ultrasonic Processing

Equipment: Ultrasonic cleaner, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5. Note: Ultrasonic cleaning may cause additional damage to implants that have surface damage.

- Prepare a solution using warm tap water and detergent or cleaner. Follow the enzymatic cleaner or detergent manufacturer's recommendations for use paying close attention to the correct exposure time, temperature, water quality, and concentration.
- Clean implants ultrasonically for a minimum of 15 minutes.
- Rinse implant thoroughly with DI or purified water.

Use DI or purified water for final rinse.

- Dry the implant using a clean soft, lint-free cloth to avoid scratching the surface.

Mechanical Processing

Equipment: Washer/disinfector, neutral enzymatic cleaner or neutral detergent with a pH ≤ 8.5.

CYCLE	MINIMUM TIME (MINUTES)	MINIMUM TEMPERATURE/WATER	TYPE OF DETERGENT
Pre-wash	2	Cold tap water	N/A
Enzyme wash	2	Warm tap water	Neutral enzymatic pH ≤ 8.5
Wash II	5	Warm tap water (>40°C)	Detergent with pH ≤ 8.5
Rinse	2	Warm DI or purified water (>40°C)	N/A
Dry	40	90°C	N/A

Instrument Cleaning Requirements: Acumed Instruments and Accessories must be thoroughly cleaned before reuse. Decontamination of reusable instruments or accessories should occur immediately after completion of the surgical procedure. Excess blood or debris should be wiped off to prevent it from drying onto the surface. Using an enzymatic cleaning product such as Enzol, prepare the cleaning

solution according to manufacturer's instructions.

Avoiding prolonged soaking, submerge the device in the enzymatic cleaning solution for 5 minutes. After 5 minutes, use a soft bristle brush to remove any detergent residue. Instruments with lumens also should be held under water (larger side up) ensuring the lumen is flushed with water. Place the device in an ultrasonic cleaner filled with deionized water, ensuring that there is no air trapped in any of the crevices by gently agitating the device. Sonicate the device for 5 minutes. Remove the device from the ultrasonic cleaner and rinse in water by submerging it completely. Gently wipe the device with a clean, dry cloth, then air dry.

Even surgical instruments manufactured from high grade stainless steel must be dried thoroughly to prevent rust formation, and all devices must be inspected for cleanliness of surfaces, joints, and lumens, proper function, and wear and tear prior to sterilization.

Instruments manufactured from different metals should be processed separately to avoid electrolytic

action between the different metals.

Exposed Springs, Coils, or Flexible Features

Flood the crevices with copious amounts of cleaning solution to flush out any soil. Scrub the surface with a scrub brush to remove all visible soil from the surface and crevices. Bend the flexible area and scrub the surface with a scrub brush. Rotate the part while scrubbing to ensure that all crevices are cleaned.

Anodized Aluminum

Anodized aluminum must not come in contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants or solutions containing iodine, chlorine or certain metal salts. Also, in solutions with pH values above 11, the anodization layer may dissolve.

STERILITY:

System components may be provided sterile or nonsterile.

Sterile Product: Sterile product was exposed to a minimum dose of 25.0-kGy gamma irradiation. Resterilization of a device provided sterile should be handled according to hospital protocol. Acumed does

not recommend resterilization of sterile-packaged product.

Non-Sterile Product: Unless clearly labeled as sterile and provided in an unopened sterile package provided by Acumed, all implants and instruments must be considered nonsterile, and sterilized by the hospital prior to use. Nonsterile devices have been validated using the sterilization parameters listed below, in fully-loaded trays with all parts placed appropriately.

Sterilization Methods

- Consult your equipment manufacturer's written instructions for specific sterilizer and load configuration instructions.
- Follow current AORN "Recommended Practices for Sterilization in Perioperative Practice Settings" and ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- Flash sterilization is not recommended, but if used, should only be performed according to requirements of ANSI/AAMI ST79: 2010 – Comprehensive guide to steam sterilization and sterility assurance in health

care facilities.

80-1325/80-1326; 80-1332/80-1333; 80-1337/80-1338	
These trays have been validated to the parameters below:	
Gravity Displacement Autoclave: NOT RECOMMENDED	
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270° F (132° C)
Exposure Time:	4 minutes
Dry Time:	30 minutes

STORAGE INSTRUCTIONS: Store in a cool dry place and keep away from direct sunlight. Prior to use, inspect product package for signs of tampering, or water contamination. Use oldest lots first.

SYMBOL LEGEND	
	Consult instructions for use
	Caution
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Use-by date
	Batch code
	Authorized representative in the European Community
	Manufacturer
	Date of manufacture
	Do not re-sterilize
	Do not re-use
	Upper limit of temperature

APPLICABILITY: These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries.

Nothing contained on these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way which is not authorized under the laws and regulations of the country where the reader is located.

FURTHER INFORMATION: To request further material, please see the contact information listed on this document.

Cautions: Federal Law (USA) restricts this product sale by or on the order of a physician or hospital. Professional Use Only.

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INDICATIONS: The Acumed Pelvic Bone Plate System includes plates, screws and accessories designed to provide fixation during fractures, fusions, and osteotomies for the acetabulum, sacrum, ilium, and entire pelvic ring, as well as treatment of sacroiliac joint dislocations and symphysis pubis disruptions.

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1 Shellock, F. G. Reference Manual for Magnetic Resonance Safety, Implants, and Devices: 2011 Edition. Biomedical Research Publishing Group, 2011.

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CYCLE	MINIMUM TIME (MINUTES)	MINIMUM TEMPERATURE/WATER	TYPE OF DETERGENT
Pre-wash	2	Cold tap water	N/A
Enzyme wash	2	Warm tap water	Neutral enzymatic pH ≤ 8.5
Wash	5	Warm tap water (>40°C)	Detergent with pH ≤ 8.5
Rinse	2	Warm DI or purified water (>40°C)	N/A
Dry	40	90°C	N/A

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care facilities.

80-1325/80-1326; 80-1332/80-1333; 80-1337/80-1338	
These trays have been validated to the parameters below:	
Gravity Displacement Autoclave:	NOT RECOMMENDED
Pre-Vacuum Autoclave:	
Condition:	Wrapped
Exposure Temperature:	270o F (132o C)
Exposure Time:	4 minutes
Dry Time:	30 minutes

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OBSOLETE
 Visit www.aelumed.net/ifu for the latest version.

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