

# NEURO PACK

Sterile Neuro Fixation System









*Rethinking Possibilities, Reshaping Lives*




# NEURO PACK

Sterile Neuro Fixation System

- OsteoMed offers NeuroPack™ kit combinations for cranial flap fixation.
- All items are contained in one convenient kit, making implant selection quick and easy
- NeuroPacks include OsteoMed's patent Auto-Drive self drilling screws, which eliminate the need for pre-drilling.
- Complements Existing Fast-Flap™ and Profile 0™ Systems



Standard Fast Flap		
218-0070-SP		NeuroPack, 1 Long & 2 Medium Straight plate & 6 AutoDrive 4mm Screws
218-0074-SP		NeuroPack, 3 Medium Straight plate & 6 AutoDrive 4mm Screws
218-0076-SP		NeuroPack, 2 Medium Straight plate & 4 AutoDrive 4mm Screws
218-0078-SP		NeuroPack, 1 Large Bur Hole plate & 6 AutoDrive 4mm Screws
218-0083-SP		NeuroPack, 1 Bur Hole Shunt Plate with Groove & 6 AutoDrive 5mm Screws

Low Profile Fast Flap		
218-0174-SP		NeuroPack, 3 LP Medium Straight plate & 6 AutoDrive 4mm Screws
218-0176-SP		NeuroPack, 2 LP Medium Straight plate & 4 AutoDrive 4mm Screws
218-0178-SP		NeuroPack, 1 LP Large Bur Hole plate & 6 AutoDrive 4mm Screws

Profile 0		
218-0075-SP		NeuroPack, Profile Zero, 3 straight plates & 6 AutoDrive 4mm screws
218-0077-SP		NeuroPack, Profile Zero, 2 straight plates & 4 AutoDrive 4mm screws
218-0079-SP		NeuroPack, Profile Zero, 1 Large Bur Hole Plate & 6 AutoDrive 4mm screws
218-0084-SP		NeuroPack, Profile Zero, 1 2 Hole Plate & 2 AutoDrive 4mm screws
218-0086-SP		NeuroPack, Profile Zero, 1 Square Plate & 4 AutoDrive 4mm screws
218-0087-SP		NeuroPack, Profile Zero, 1 4 Hole Long Plate & 2 AutoDrive 4mm screws
218-0088-SP		NeuroPack, Profile Zero, 1 Double Y Hole Plate & 4 AutoDrive 4mm screws
218-0089-SP		NeuroPack, Profile Zero, 1 2x3 Quad Plate & 4 AutoDrive 4mm screws
218-0098-SP		NeuroPack, Profile Zero, Bur Hole Shunt Plate & 3 AutoDrive 4mm screws



3885 Arapaho Rd.  
Addison, TX 75001, USA  
Phone: +1 972 677 4600  
www.osteomed.com