


INTERPHLEX™

FLEXIBLE STABILIZATION RODS



- Available in 2 system sizes
 - Metatarsalphalangeal Joint
 - Intrapalangeal Joint
- Provides joint space stability and reduces toe migration
- Maintains flexibility of the joint
- Minimizes inflammatory response



Rethinking Possibilities, Reshaping Lives

INTERPHLEX™

FLEXIBLE STABILIZATION ROD

The **OsteoMed InterPhlex System** is an alternative to arthroplasty or arthrodesis of the lesser joints. The double stemmed implant is designed to maintain flexibility, length and shape of the toe. Made from medical grade silicone the rods have a bulb that provides stability of the joint space while the stiffness of the rod addresses the concern for toe migration.



INTER-PHALANGEAL

360-2200 Inter-Phalangeal Sterile Surgical Kit, Drills & Trials



METATARSAL PHALANGEAL

360-2800 Metatarsal Phalangeal Sterile Surgical Kit, Drills & Trials

OsteoMed Products



ExtremiLOCK Ankle Plating System



ExtremiLOCK Ankle Plating System



ExtremiLOCK Foot Plating System



ExtremiFix Small Cannulated Screws



ExtremiFix LDA Screws



Large Cannulated Screws



M3-X



EnCompass



ReFlexion



Hemi



ExtremiFuse



Talar-Fit



KobyGard



Inion



OsteoVation EX



OsteoVation QWIK



OsteoVation B-TCP



Allograft Offering