


OSTEOMED

A TRADITION OF QUALITY AND INNOVATION

OSTEOFORM™

FORMABLE CRANIOFACIAL MESH



STRENGTH AND PASSIVE CONTOURABILITY THROUGH ENGINEERING

- PRE-FORMED MESH ALLOWS FOR EASY ADAPTATION TO ANATOMICAL CURVATURE OF CRANIUM
- APPLICABLE TO VARYING ANATOMICAL CURVES WITH MINOR CONTOURING OF APPLIANCE
- RAPID MESH APPLICATION WHEN USED WITH OSTEOMED'S PATENTED AUTO-DRIVE™ SELF-DRILLING SCREWS
- RIGIDITY AND STRENGTH PROVIDE FOR BETTER PROTECTION OF UNDERLYING STRUCTURES
- USE AS A STAND ALONE TITANIUM APPLIANCE OR AS SCAFFOLDING FOR CRANIOPLASTY


PRODUCTS THAT IMPROVE PEOPLE'S LIVES

Patent Number: 5,925,048


OSTEOFORM™

FORMABLE CRANIOFACIAL MESH

Refer to *Product Insert* for specific indications and *Instructions For Use*.



Unique Mesh Designs



Anatomically
Preformed

<u>PART NUMBER</u>	<u>DESCRIPTION</u>
218-2001	40mm X 40mm Malleable OsteoForm Mesh Plate
218-2002	40mm X 40mm Rigid OsteoForm Mesh Plate
218-2003	90mm X 90mm Malleable OsteoForm Mesh Plate
218-2004	90mm X 90mm Rigid OsteoForm Mesh Plate
218-2005	165mm X 95mm Pre-Formed OsteoForm Mesh Plate
218-2006	90mm X 90mm Pre-Formed OsteoForm Mesh Plate
218-2007	165mm x 95mm Rigid OsteoForm Mesh Plate
218-2008	165mm x 95mm Pre-Formed OsteoForm Mesh, Extra-Rigid

218-2001-SP	40mm x 40mm Malleable OsteoForm Mesh, Sterile
218-2002-SP	40mm x 40mm Rigid OsteoForm Mesh, Sterile
218-2003-SP	90mm x 90mm Malleable OsteoForm Mesh, Sterile
218-2004-SP	90mm x 90mm Rigid OsteoForm Mesh, Sterile
218-2005-SP	165mm x 95mm Pre-Formed OsteoForm Mesh, Sterile
218-2006-SP	90mm x 90mm Pre-Formed OsteoForm Mesh, Sterile



3885 Arapaho Road
Addison, Texas 75001
Customer Service: 800.456.7779
Main: 972.677.4600
Fax: 800.390.2620

www.osteomedcorp.com

P/N 030-1243 Rev. D